

Internal SD Controller Backup

Materials

- USB thumb drive

Protocol

1. Locate the back access panel on the back of the KinExA® instrument (**Figure 1**). If there is no USB port on your panel, contact Sapidyne for assistance.
2. With the instrument turned on, plug a USB thumb drive into the USB port indicated in **Figure 1**.
3. The busy (orange) LED light on the front panel will illuminate once the instrument recognizes the USB thumb drive.
 - If the Error (red) LED light illuminates, then the USB thumb drive does not have the necessary capacity to store the backup. You may need up to 2 gigabytes of free space on the thumb drive.
4. The busy (orange) LED light will remain illuminated while the files are being uploaded.
5. Once the files are uploaded, the busy (orange) LED light will blink indicating that the files have been successfully uploaded. It is now safe to remove the USB thumb drive from the USB port.




Figure 1. Back access panel.